



FREE Flexi-Block.
With each Techne® Dri-Block® heater.
Promotion until 30.06.2015.

The new Techne® Flexi-block accessory is ideal for heating samples stored in 0.2ml, 0.5ml, 2.0ml and 5.0ml micro tubes, as well as any diameter test-tube, vial, cuvette or micro plate. Manufactured from high density ceramic, the autoclaveable beads exhibit high mechanical strength and are an excellent thermal transfer medium.



FREE Flexi-Block
worth £190

FREE FLEXI2 block insert.

Please use code TH115OFFER on your purchase order



FREE Flexi-Block
worth £285

FREE FLEXI3 block insert.

Please use code TH115OFFER on your purchase order



www.techne.com



FREE Flexi-Block.

With each Techne® Dri-Block® heater.

Promotion until 30.06.2015.

Dri-Block® heater Key Features

- Integrated count-down timer, 1m to 99h59m with alarm
- New fascia with brighter 5 character LED Display
- Small, light compact footprint; can be stacked in storage
- Fast heating: ambient to 100°C in 10 minutes
- 3 year warranty

Flexi-Block® Key Features

- Heat any sample
- Block made from aluminium
- Autoclaveable ceramic beads included
- Suitable for all Techne® 2 and 3 Dri-block® units
- Manufactured in UK



FREE FLEXI2 Flexi-Block

Part number	Description	Price
DB100/2	Dri-Block® heater. Ambient 5°C to 100°C, with FREE FLEXI2	£610.00
FLEXI2	Flexi-block, double-width. Heats tubes, vials and unskirted plates. Beads included.	FREE

Please use code TH115OFFER on your purchase order.

FREE FLEXI3 Flexi-Block

Part number	Description	Price
DB100/3	Dri-Block® heater. Ambient 5°C to 100°C, with FREE FLEXI3	£703.00
FLEXI3	Flexi-block, triple-width. Heats tubes, vials and skirted plates. Beads included.	FREE

Please use code TH115OFFER on your purchase order.

Terms and Conditions

Valid until 30.06.2015. All prices subject to standard terms and conditions of sale. Special offer code MUST be used at time of ordering to receive the offer.



info@techne.com

www.techne.com

