

Elite Thermal Cyclers. **Sale until 31/07/14.**
Available from your **Techne®** Distributor.



Prime Elite special price
£3200 / €3840 / \$5120

Please use code ELITE1 on your purchase order

Prime Satellite special price
£2500 / €3020 / \$4030

Please use code ELITE1 on your purchase order

Grow your network Techne Prime Elite Advanced Cycler

This premium thermal cycler is equipped as standard with a large colour touchscreen, large 30°C temperature gradient, fast 5°C/sec ramp rates, an innovative front-loading PCR plate drawer, TERSTM thermal energy recovery system and a unique space-saving stackable design. Prime Elite and Satellites combine to form a multi-block network. Satellites are controlled directly from a Prime Elite or from a PC with the supplied software.



Prime Elite + 1 Satellite
£5500 / €6600 / \$8800

Please use code ELITE2 on your purchase order

Prime Elite + 2 Satellite
£7700 / €9240 / \$12320

Please use code ELITE3 on your purchase order

Prime Elite + 3 Satellite
£9700 / €11640 / \$15520

Please use code ELITE4 on your purchase order



www.thebestthermalcycler.com



Choose your capacity.


1 Prime Elite	£3200 / €3840 / \$5120	Promo Code ELITE1
1 Prime Elite Satellite*	£2500 / €3020 / \$4030	Promo Code ELITE1
1 Prime Elite + 1 Satellite	£5500 / €6600 / \$8800	Promo Code ELITE2
1 Prime Elite + 2 Satellites	£7700 / €9240 / \$12320	Promo Code ELITE3
1 Prime Elite + 3 Satellites	£9700 / €11640 / \$15520	Promo Code ELITE4

* Satellites are controlled directly from a Prime Elite or from a PC with supplied software.

Ordering information

Product Code	Description	Please note
ELITE/02	Prime Elite gradient thermal cycler, 96 x 0.2ml	For dual, triple or quad systems
ELITESAT/02	Prime Elite satellite thermal cycler, 96 x 0.2ml	Please order multiples satellites

Technical Specification

	Prime Elite		Prime Elite Satellite
Sample capacity: 0.2ml	96		96
Maximum heating rate	5.0°C/sec		5.0°C/sec
Block temperature range	4°C to 100°C		4°C to 100°C
Block uniformity at 55°C	±0.3°C		±0.3°C
Temperature accuracy at 55°C	±0.2°C		±0.2°C
Gradient calculator	Integrated		Via Prime Elite or PC
Type	8 peltier elements/block		8 peltier elements/block
Selectable heated lid temp	35°C to 115°C or off		35°C to 115°C or off
Program interface	5.7" VGA touchscreen		Via Prime Elite or PC
Number of programs	1000		Via Prime Elite or PC
Data transfer	USB port		USB port
Auto re-start on power failure	Yes		Yes
Dimensions (w x d x h), mm	276 x 375 x 312mm		276 x 375 x 312mm
Weight, kg	14		14
Voltage	100-230V, 50-60Hz		100-230V, 50-60Hz
Power	1000W		1000W
Country of manufacture	UK		UK
Warranty	4 years		4 years



Part of the Techne® PRIME Range
Find out more information via:

info@techne.com

www.techne.com